

Set per Estrazione

Dove nasce un sorriso...
Where smile rises...



- ▶ **Set Strumenti per Estrazione e per la Chirurgia Orale**
- ▶ *Instruments SET for Extraction and Oral Surgery*



S0100-1

Contenuto/Contents: 0100-1, 0100-17, 0100-18, 0100-79, 0100-67A, 0100-51A, 0100-13, 0100-22.



S0100-2

Contenuto/Contents: 0100-1, 0100-7, 0100-17, 0100-18, 0100-4, 0100-13, 0100-22, 0100-86C.



S0100-3

Contenuto/Content: 0100-67A, 0100-51A, 0100-33, 0100-79, 0200-1, 0201-2, 0201-3, 0202-320, 0202-321.



S0100-4

Contenuto/Contents: 0100-67A, 0100-51A, 0100-33, 0100-79, 0200-1, 0201-2, 0201-3, 0435-14TC, 0425-16TC.

8002

Steriloil - Contenuto/Contents: 40 ml



Massima protezione degli strumenti in autoclave. Penetra negli snodi e lubrifica. Protegge e lubrifica gli strumenti chirurgici e dentali prima e dopo la sterilizzazione. Modalità d'uso: oliare prima e dopo la sterilizzazione le punte, gli snodi e le chiusure degli strumenti.

Best instruments' autoclaving protection. Penetrates and lubricates all instruments' joints. Protects all dental and surgical hand instruments before and after sterilization. Instructions for use: oil all tips, joints and fastening of instruments before and after sterilization.

8010

Silicone - Contenuto/Contents: 40 ml



Lubrificante speciale per strumenti dentali e chirurgici. Per lubrificare gli snodi di forbici, pinze per estrazione, ossivore, emostatiche, porta aghi ecc. prima e dopo la sterilizzazione. Non prosciuga durante la sterilizzazione. Protegge dalla corrosione e da altri depositi e cura il buon funzionamento degli strumenti.

Specific lubricant for dental and surgical instruments. The lubricant is designed to lubricate joints of all dental instruments (scissors, haemostatic and extracting forceps, bone rongeurs, needle holders etc.). Before and after sterilization it does not dry out during sterilization. Provides protection from corrosion and other deposits and ensures proper operation of the instruments.